
#Easy Handling & Moldable **#Chair-Time Saving** #Space-Making Stability & Volume Maintenance









© Purgo Biologics Solution














Manufacturer Purgo Biologics Inc. 812, 27, Dunchon-daero 457beon-gil, Jungwon-gu, Seongnam-si, Gyeonggi-do, 13219, Korea Tel. +82 70 4827 5352 | E-mail. overseas@purgobio.com | www.purgobio.com



#### Contents

#### I. Article


- 1. Augmentation Stability of Guided Bone Regeneration for Peri-Implant
  Dehiscence Defects with L-shaped Porcine-Derived Block Bone Substitute
- 2. Use of Porcine-derived Block Bone Substitutes for Guided Bone Regeneration in the Peri-implant Dehiscence Defects of the Mandibular Anterior Region
- 3. Profilometric, volumetric, and esthetic analysis of guided bone regeneration with L-shaped collagenated bone substitute and connective tissue graft in the maxillary esthetic zone: A case series with 1-year observational study

### — II. GBR with THE Graft Collagen in Peri-implant Dehiscence Defect

| 1. Case 01    | 80 |
|---------------|----|
| 2. Case 02    | 10 |
| 3. Conclusion | 12 |

#### III. Clinical Cases

| 1. Extraction sockets                                  | 14 |
|--|----|
| 2. Horizontal bone defect                              | 24 |
| 3. Vertical and horizontal bone defect                 | 28 |
| 4. Sinus pneumatization and reduced residual bone loss | 32 |
| 5. Peri-implantitis                                    | 34 |



#### **Article**

#### 01

## Augmentation Stability of Guided Bone Regeneration for Peri-Implant Dehiscence Defects with L-shaped Porcine-Derived Block Bone Substitute

#### Jae-Hong Lee | Eun-Hee Jung | Seong-Nyum Jeong

• Department of Periodontology, Daejeon Dental Hospital, Institute of Wonkwang Dental Research, Wonkwang University College of Dentistry, Daejeon, Korea

Purpose: Block bone substitutes have better augmentation stability for guided bone regeneration (GBR) than particulate bone substitutes. This study sought to determine whether GBR with an L-shaped porcine block bone (DPBM-C; THE Graft Collagen) differs from GBR with an L-shaped bovine block bone (DBBM-C) based on clinical, radiographic, and volumetric outcomes for peri-implant dehiscence defects. A total of 42 peri-implant defects were grafted with 20 L-shaped DPBM-C and 22 DBBM- C groups.

**Method:** The horizontal and vertical thicknesses of the augmented hard tissue were measured using sagittal cone-beam computed tomography, and the volumetric tissue change was evaluated by stereolithography image superimposition. Postoperative discomfort, early wound healing outcomes, and implant stability were also assessed.

Results: Among the clinical (subjective pain and swelling, wound dehiscence, membrane exposure, and periotest values),radiographic(changes in horizontal and vertical hard tissue thickness), and volumetric parameters of the L-shaped DPBM-C and DBBM-C groups during the healing period, only the periotest values showed a statistically significant difference (0.67 1.19, p = 0.042).

Conclusion: Within the limitations of this study, an L-shaped DPBM-C is not inferior to an L-shaped DBBM-C based on their clinical, radiographic, and volumetric outcomes for GBR of peri-implant dehiscence defects.



LEE, Jae-Hong; JUNG, Eun-Hee; JEONG, Seong-Nyum. Augmentation stability of guided bone regeneration for peri-implant dehiscence defects with L-shaped porcine-derived block bone substitute. Materials, 2021, 14.21: 6580.

#### 02

# Use of Porcine-derived Block Bone Substitutes for Guided Bone Regeneration in the Peri-implant Dehiscence Defects of the Mandibular Anterior Region

#### Jae-Hong Lee | Seong-Nyum Jeong

• Department of Periodontology, Daejeon Dental Hospital, Institute of Wonkwang Dental Research, Wonkwang University College of Dentistry, Daejeon, Korea

**Purpose:** Various bone graft materials are currently used for guided bone regeneration (GBR), but data are lacking on the superiority of one biomaterial relative to the other. In addition, achieving successful bone augmentation remains challenging. When used in GBR, block bone substitutes show better augmentation stability of the bone than particulated bone substitutes. Recently developed soft-type block bone substitutes are widely used in clinical practice.

Method: : This study enrolled 12 patients who underwent implant placement with GBR in the mandibular anterior region between September 2016 and May 2021 in the Department of Periodontology, Daejeon Dental Hospital, Wonkwang University. GBR was performed using demineralized porcine bone mineral with 10% collagen (DPBM-C; THE Graft® Collagen) and a resorbable collagen membrane. The surgical site was sutured with 4-0 e-PTFE (Biotex®) and 5-0 nylon (Monosof®) by using modified horizontal mattress and interrupted suture methods. Stitch-out was performed two weeks after the operation, and reentry surgery was performed within four to six months after the implant surgery.

Results: There was a significant decrease in bone thickness after five months of implant surgery with GBR compared with the decrease immediately after surgery. However, bone gains of more than 1.5 mm were observed at all measurement sites compared with those at baseline.

Conclusion: GBR performed using demineralized porcine bone mineral with 10% collagen in combination with a resorbable collagen membrane is an effective and favorable treatment modality for peri-implant dehiscence defects in terms of clinical and radiographic outcomes during the healing period.





#### 03

Profilometric, volumetric, and esthetic analysis of guided bone regeneration with L-shaped collagenated bone substitute and connective tissue graft in the maxillary esthetic zone: A case series with 1-year observational study

#### Jae-Hong Lee | Eun-Hee Jung | Seong-Nyum Jeong

• Department of Periodontology, Daejeon Dental Hospital, Institute of Wonkwang Dental Research, Wonkwang University College of Dentistry, Daejeon, Korea

Purpose: The aim of this study was to evaluate 1-year stability and maintenance of peri-implant soft and hard tissues after guided bone regeneration (GBR) with L- shaped collagenated bone substitute (THE Graft Collagen) and subepithelial connective tissue graft (CTG) in the maxillary anterior region using profilometric, volumetric, and esthetic analyses.

Method: Fourteen peri-implant defects were grafted with L-shaped collagenated bone substitute, and 5 months after implant placement with GBR, reentry surgery in combination with CTG was performed in all participants. CBCT scans and STL files were acquired at baseline (after implant surgery, T1), reentry surgery (T2), and 1-year follow-up (T3). The profilometric and volumetric changes of the peri-implant tissues were measured, and the pink esthetic score (PES) was assessed at T3.

Results: One year after GBR and CTG at the buccal aspect of the maxillary esthetic zone, the mean thickness of the hard tissue (HT) decreased (HT0:  $0.87 \pm 0.67$  mm, HT1:  $0.74 \pm 0.75$  mm, HT2:  $0.92 \pm 0.48$  mm, 45-HT:  $0.87 \pm 0.73$  mm) and the corresponding thickness of the soft tissue (ST) increased (ST0:  $0.96 \pm 1.06$  mm, ST1:  $0.85 \pm 0.95$  mm, ST2:  $0.38 \pm 0.82$  mm, 45-ST:  $0.12 \pm 0.62$  mm), and as a result, there was no statistically significant difference in the total tissue thickness between T1 and T3 (p < 0.05). The mean volumetric changes of the peri-implant tissues increased after 1-year of implant surgery (T1–T2:  $1.52 \pm 0.83$  mm, T2–T3:  $0.88 \pm 1.04$  mm, T1–T3:  $0.64 \pm 0.90$  mm), and a statistically significant difference was shown in all compared time periods (p < 0.05). The mean PES score was  $8.07 \pm 1.54$  at T3 (range, 6–10).

Conclusion: Within the limitations of this 1-year follow-up study, GBR with an L- shaped collagenated bone substitute and subepithelial CTG in the maxillary esthetic zone was beneficial for stable and maintainable peri-implant hard and soft tissues.



LEE, Jae-Hong; JUNG, Eun-Hee; JEONG, Seong-Nyum. Profilometric, volumetric, and esthetic analysis of guided bone regeneration with L-shaped collagenated bone substitute and connective tissue graft in the maxillary esthetic zone: A case series with 1-year observational study. Clinical Implant Dentistry and Related Research, 2022, 24.5: 655-663.



#### **DEHISCENCE DEFECT** | BONE AUGMENTATION

#### Case01


## Mandibular peri-implant dehiscence defect. GBR with simultaneous implant placement

#### Jung-Seok Lee

• Department of Periodontology, Research Institute for Periodontal Regeneration, Yonsei University College of Dentistry


#### **Case Summary** Nationality Korean Late-forties Grafting Simultaneous approach: Implant placement (two-stage) and GBR with THE Treatment Plan Graft Collagen and resorbable collagen membrane. Category Products Method Description of the method Bone graft THE Graft Collagen Single used Materials A buccal dehiscence defect following implant placement was augmented using The Graft Collagen and a resorbable collagen membrane. Resorbable collagen Membrane membrane Approach Category Indication Surgical Procedures ☐ Alveolar ridge ■ Implantology □ Extraction sockets □ One-stage □ Periodontology ■ Dehiscence □ Intra-socket ■ Two-stage ☐ Fenestration ■ Bone augmentation ☐ Immediate placement/Immediate loading ■ Horizontal bone loss ☐ Ridge Split Methods ■ Simultaneous approach □ Vertical bone loss ☐ Staged/Delayed approach □ Sinus pneumatization □ Peri-implantitis □ Furcation

#### **Before**




Observation of damaged labial alveolar bone.

#### **After**



1 year post-GBR. Observation of corticalization similar to the adjacent alveolar bone on radiographs.



1. Pre-operation

2. Pre-operation, CBCT findings. (a) Axial view; (b) Coronal view. Observation of damaged labial


3. Intra-first operation. Incision and flap elevation.



 Intra-first operation. Implant placement and THE Graft Collagen were applied in an L-shape to the buccal and occlusal surface.

5. Intra-first operation. Covered with a resorbable

6. Intra-first operation. Sutured




7. Post-first operation, CBCT findings. (a) Axial view; 8. Pre-second operation, 4 months post-GBR. (b) Coronal view. Observation of increased volume






9. Pre-second operation, 4 months post-GBR.



10. Post-second operation, 4 months post-GBR. (a) Axial view; (b) Coronal view. Natural bone remodeling occurs according to the bone envelope.



11. Final prosthesis, 1 year post-GBR.



12. Final prosthesis, 1 year post-GBR. (a) Axial view; (b) Coronal view. Observation of corticalization similar to the adjacent alveolar bone on radiographs.

#### **DEHISCENCE DEFECT** | BONE AUGMENTATION

## Case02


## Mandibular peri-implant dehiscence defect. GBR with simultaneous implant placement

#### Dong-Woon Lee

• Department of Periodontology, Dental Hospital, Wonkwang University College of Dentistry


| Nationality    | Korean  | Age                          | Mid-Seventies                 |  |                               |  |
|----------------|---|------------------------------|-------------------------------|--|-------------------------------|--|
| Treatment Plan | Simultaneous approach: I<br>Graft Collagen and resort |                              |                               | Grafting<br>Area   | 6                             |  |
|                | Category  | Products                     | Method                        | D  | escription of the method      |  |
| Materials      | Bone graft  | THE Graft Collagen           | Single used                   | A buccal dehiscence defect following implacement was augmented using |                               |  |
|                | Membrane  | Resorbable collagen membrane |                               | THE Graft Collagen and a resorbable collag membrane.                 |                               |  |
|                | Category  | Indication                   | Approach                      |  | Surgical Procedures           |  |
|                | ■ Implantology  | ☐ Extraction sockets         | ☐ Alveolar ridge preservation | □ One-stag   | е                             |  |
|                | □ Periodontology                                      | ■ Dehiscence                 | □ Intra-socket                | ■ Two-stage  | е                             |  |
|                |   | □ Fenestration               | ■ Bone augmentation           | □ Immediate  | e placement/Immediate loading |  |
| Methods        |   | ■ Horizontal bone loss       | □ Ridge Split                 | ■ Simultane  | eous approach                 |  |
|                |   | □ Vertical bone loss         | □ Lateral                     | □ Staged/D   | elayed approach               |  |
|                |   | ☐ Sinus pneumatization       | □ Crestal                     |  |                               |  |
|                |   | □ Peri-implantitis           |                               |  |                               |  |
|                |   | □ Furcation                  |                               |  |                               |  |

#### **Before**




Observation of damaged buccal alveolar bone.

#### **After**



4 months post-GBR. Natural bone remodeling occurs according to the bone envelope.



2. Pre-operation, CBCT findings. (a) Axial view; (b) Coronal view. Observation of damaged labial alveolar bone.

3. Intra-first operation. Incision and flap elevation.




4. Intra-first operation. Implant placement.

 Intra-first operation. THE Graft Collagen is applied
in an L-shape to the buccal and occlusal surfaces
 Intra-first operation. Sutured. and covered with a collagen membrane.



7. Post-first operation, CBCT findings. (a) Axial view; 8. Pre-second operation, 4 months post-GBR. (b) Coronal view. Observation of increased buccal bone volume following GBR.

9. Intra-second operation, 4 months post-GBR.



10. Post-second operation, 4 months post-GBR.

11. Post-second operation, 4 months post-GBR.
(a) Axial view; (b) Coronal view. Natural bone remodeling occurs according to the bone envelope.

12. Final prosthesis, 9 months post-GBR.

#### Conclusion

## Regeneration of peri-implant dehiscence defect resolution by GBR with simultaneous implant placement


Prof. Jung-Seok Lee<sup>1</sup> | Prof. Dong-Woon Lee<sup>2</sup>

- 1 Department of Periodontology, Research Institute for Periodontal Regeneration, Yonsei University College of Dentistry
- 2 Department of Periodontology, Dental Hospital, Wonkwang University College of Dentistry

A block-type collagenated bone substitute cannot be categorized as a conventional 'block-bone' that is widely used for regeneration of an extensive alveolar bone defect. Even though it has the block morphology and is able to provide cohesion of particulate bone substitutes at the time of surgery, the original shape and volume may be lost at an early stage by the rapid resorption of collagen binder within the material. In order to enhance the space-maintaining property of the block, THE Graft Collagen has been manufactured with a binder of physically crosslinked collagen matrix, which is able to effectively maintain the augmented volume by slowing down the resorption rate of the collagen binder.

In the clinical cases in this case-book, THE Graft Collagen had been grafted beyond the alveolar bony envelope, considering the possible resorption of the collagen binder. Four months later, we observed restoration of the alveolar ridge with bone regeneration as well as substantial shrinkage of the augmented volume. Thus, alveolar ridge regeneration could be achieved up to the bony envelope extending from the outer shape of the adjacent alveolar ridge. Based on the reference line of bony envelope, the conventional guided bone regeneration (GBR) and the simplified GBR with THE Graft Collagen could produce similar bone regeneration. In most cases, the reconstruction of alveolar ridge up to the same level as the adjacent tissues is clinically sufficient to surround and support the dental implant. Therefore, in these sites, THE Graft Collagen could be a good option for successful regeneration of alveolar ridge. In addition, the collagen binder enhances manageability of the grafted materials and maintains the bone substitutes within the defect site. THE Graft Collagen could be a successful option for bone grafting technique due to its simplicity to apply in the defect and capacity to maintain its stability.

This work was supported by the Korea Medical Device Development Fund grant funded by the Korea government (the Ministry of Science and ICT, the Ministry of Trade, Industry and Energy, the Ministry of Health & Welfare, the Ministry of Food and Drug Safety) (Project number: ,KMDF\_PR\_20200901)



# Maxillary right and left central incisor implants with intra-socket bone graft and ridge preservation of maxillary right and left lateral incisors

#### Yong-Seok Cho


• The 22<sup>nd</sup> Century Seoul Dental Hospital



#### **Case Summary**

| Nationality     | Korean  | Age  | Late Fifties                |  |  |                  |         |
|-----------------|---|--|-----------------------------|--|--|------------------|---------|
| Chief Complaint |   | uxillary right and left lateral in of maxillary right and left c |                             |  |  |                  |         |
| Treatment Plan  | Extraction of the upper lateral and central incisors, bilaterally.     Immediate implant on placement on the left and right central incisor areas.     Alveolar ridge preservation using a collagenated bone block. |  |                             |  | 1. Extraction of the upper lateral and central incisors, bilaterally.     2. Immediate implant on placement on the left and right central incisor areas. | Grafting<br>Area | 2 1 1 2 |
|                 | Category  | Products   | Method                      | De   | scription of the method  |                  |         |
| Materials       | Bone graft  | THE Graft Collagen   | Single used                 | Alveolar ridge preservation and intra-socke filling. |  |                  |         |
|                 | Category  | Indication   | Approach                    |  | Surgical Procedures  |                  |         |
|                 | ■ Implantology  | ■ Extraction sockets   | Alveolar ridge preservation | ■ One-stage  | 9  |                  |         |
|                 | □ Periodontology  | □ Dehiscence   | ■ Intra-socket              | ☐ Two-stage  | •  |                  |         |
|                 |   | ☐ Fenestration   | ☐ Bone augmentation         | ■ Immediate  | e placement/Immediate loading  |                  |         |
| Methods         |   | ☐ Horizontal bone loss   | ☐ Ridge Split               | ■ Simultane  | ous approach   |                  |         |
|                 |   | ☐ Vertical bone loss   | □ Lateral                   | □ Staged/De  | elayed approach  |                  |         |
|                 |   | ☐ Sinus pneumatization   | □ Crestal                   |  |  |                  |         |
|                 |   | □ Peri-implantitis   |                             |  |  |                  |         |
|                 |   | □ Furcation  |                             |  |  |                  |         |

#### **Before**



The patient visited for the mobility of the maxillary anterior bridge.

#### **After**



The periapical X-ray showed well-integrated implants.

#### Conclusion

THE Graft Collagen made intra-socket bone grafting and ridge preservation simple and effective.



 Pre-operation. The patient visited for the mobility of the maxillary anterior bridge.

2. Pre-operation

 Intra-operation. Observed insufficient residual alveolar bone width and height.



 Intra-operation. THE Graft Collagen was used to preserve the ridges of the maxillary right and left lateral incipera.

Intra-operation. THE Graft Collagen was used to perform an intra-socket bone graft at the maxillary right and left central incisors labial gap.

6. Intra-operation. Flap approximation with 4-0 nylon.



7. Post-operation.

8. Final prosthesis, 2 months and 3 weeks post-operation.

9. Final prosthesis, 2 months and 3 weeks post-operation.



10. Final prosthesis, 2 months and 3 weeks post-operation.

11. Follow-up, 10 months post-operation.

12. Follow-up, 10 months post-operation.

## Implant with intra-socket bone graft on the maxillary left canine

#### Yong-Seok Cho

• The 22<sup>nd</sup> Century Seoul Dental Hospital



#### **Case Summary**

| Nationality     | Korean  | Age                    | Late Eighties                 |  |                               |
|-----------------|---|------------------------|-------------------------------|--|-------------------------------|
| Chief Complaint | Severe mobility of a bridge spanning from the upper left canine to the second premolar.   |                        |                               |  |                               |
| Treatment Plan  | Extraction of the maxillary left canine and second premolar.     Immediate implant placement on the maxillary left canine and second premolar.     Intra-socket bone graft with THE Graft Collagen. |                        |                               | Grafting<br>Area                               | 3                             |
|                 | Category  | Products               | Method                        | De   | escription of the method      |
| Materials       | Bone graft  | THE Graft Collagen     | Single used                   | Intra-socket bone graft and ridge preservation |                               |
|                 | Category  | Indication             | Approach                      |  | Surgical Procedures           |
|                 | ■ Implantology  | ■ Extraction sockets   | ☐ Alveolar ridge preservation | ■ One-stag                                     | е                             |
|                 | □ Periodontology  | □ Dehiscence           | ■ Intra-socket                | ☐ Two-stage                                    | е                             |
|                 |   | ☐ Fenestration         | ☐ Bone augmentation           | ■ Immediate                                    | e placement/Immediate loading |
| Methods         |   | ☐ Horizontal bone loss | □ Ridge Split                 | ■ Simultane                                    | eous approach                 |
|                 |   | ☐ Vertical bone loss   | □ Lateral                     | □ Staged/D                                     | elayed approach               |
|                 |   | ☐ Sinus pneumatization | □ Crestal                     |  |                               |
|                 |   | □ Peri-implantitis     |                               |  |                               |
|                 |   | □ Furcation            |                               |  |                               |

#### **Before**



The patient visited for discomfort on the maxillary left canine to the second premolar bridge.

**After** 

The periapical X-ray showed well-integrated implants with stable marginal bone.

#### Conclusion

Intra-socket bone graft using THE Graft Collagen was effective in preserving the alveolar bone contour.

#### **Case Summary**




- Pre-operation. The patient visited for the severe mobility of the maxillary left to the second premolar bridge.
- 2. Pre-operation.

3. Intra-operation. Maxillary left canine and second premolar were extracted. Sequential osteotomy was performed at the palatal slope of the extraction sockets.



- 4. Intra-operation. Implants were placed on the palatal slope. A > 3.0 mm gap was found on the labial
   5. Intra-operation. Bone grafting was performed to fill the gaps using THE Graft Collagen. aspects of the placed implants.
- Intra-operation. The healing abutments were connected.



7. Post-operation.

8. Post-operation.

9. Follow-up, 2 months post-operation.



- 10. Final prosthesis, 3 months and 2 weeks post-operation.
- 11. Final prosthesis, 3 months and 2 weeks post-operation.
- 12. Follow-up, 6 months and 3 weeks post-operation.

## Implant with intra-socket bone graft on the maxillary right lateral incisor and canine

#### Yong-Seok Cho

• The 22<sup>nd</sup> Century Seoul Dental Hospital



#### **Case Summary**

| Nationality     | Korean   | Age                          | Early fifties                 |                          |                               |  |
|-----------------|--|------------------------------|-------------------------------|--------------------------|-------------------------------|--|
| Chief Complaint | Mobility and discomfort of   | on the maxillary right later |                               |                          |                               |  |
| Treatment Plan  | Extraction of the maxillary right lateral incisor and canine.     Implant on the maxillary right lateral incisor and canine.     Intra-socket bone graft with a collagenated bone. |                              |                               | Grafting<br>Area         | 3 2                           |  |
|                 | Category   | Products                     | Method                        | De                       | escription of the method      |  |
| Materials       | Bone graft   | THE Graft Collagen           | Single used                   | Intra-socket bone graft. |                               |  |
|                 | Category   | Indication                   | Approach                      |                          | Surgical Procedures           |  |
|                 | ■ Implantology   | ■ Extraction sockets         | ☐ Alveolar ridge preservation | ■ One-stag               | е                             |  |
|                 | □ Periodontology   | □ Dehiscence                 | ■ Intra-socket                | □ Two-stage              | е                             |  |
|                 |  | □ Fenestration               | ☐ Bone augmentation           | ■ Immediate              | e placement/Immediate loading |  |
| Methods         |  | ☐ Horizontal bone loss       | ☐ Ridge Split                 | ■ Simultane              | eous approach                 |  |
|                 |  | ☐ Vertical bone loss         | □ Lateral                     | □ Staged/D               | elayed approach               |  |
|                 |  | ☐ Sinus pneumatization       | □ Crestal                     |                          |                               |  |
|                 |  | □ Peri-implantitis           |                               |                          |                               |  |
|                 |  | □ Furcation                  |                               |                          |                               |  |

#### Before




The patient visited for discomfort on the maxillary right lateral incisor and canine.

**After** 

The periapical X-ray showed well-integrated implants with stable marginal bone.

#### Conclusion

THE Graft Collagen made intra-socket bone grafting simple and effective.



 Pre-operation. The patient visited for mobility and discomfort on the maxillary right lateral incisor and canine.

2. Pre-operation


**3. Intra-operation.** Extraction of the maxillary right lateral incisor and canine.



**4. Intra-operation.** Implants were placed on the palatal slope.

 Intra-operation. THE Graft Collagen was used to perform an intra-socket bone graft at the labial gaps within the extraction sockets of maxillary right lateral incisor and canine.


**6.Intra-operation.** The healing abutments were connected. Labially exposed graft materials were left without



7. Post-operation.

8. Post-operation.

9. Final prosthesis, 4 months post-operation.



10. Final prosthesis, 4 months post-operation.

11. Final prosthesis, 4 months post-operation.

12. Follow-up, 5 months post-operation.

## Immediate implant placement with gap filling at the extraction sockets on the maxillary left central incisor and alveolar ridge preservation at the maxillary left lateral incisor

#### Yong-Seok Cho


• The 22<sup>nd</sup> Century Seoul Dental Hospital



#### **Case Summary**

| Nationality     | Korean  | Age   | Late forties  |   |
|-----------------|---|---|---|---|
| Chief Complaint | Crown fracture of the max   | xillary left lateral incisor.   |   |   |
| Treatment Plan  | Extraction of the maxillary left central and lateral incisors.     Immediate implant placement on the maxillary left central incisor.     Intra-socket bone graft and ridge preservation with a collagenated bone substitute. |   |   | Grafting 1 2  |
|                 | Category  | Products  | Method  | Description of the method   |
| Materials       | Bone graft  | THE Graft Collagen  | Single used   | Intra-socket bone graft and ridge preservation.                             |
|                 | Category  | Indication  | Approach  | Surgical Procedures   |
|                 |   |   |   |   |
|                 | ■ Implantology  | ■ Extraction sockets  | Alveolar ridge preservation                                   | ■ One-stage   |
|                 | ■ Implantology  □ Periodontology  | ■ Extraction sockets  □ Dehiscence                                      | ■ Alveolar ridge preservation ■ Intra-socket                  | ■ One-stage  □ Two-stage  |
|                 |   |   |   |   |
| Methods         |   | □ Dehiscence  | ■ Intra-socket  | □ Two-stage   |
| Methods         |   | □ Dehiscence □ Fenestration   | ■ Intra-socket  □ Bone augmentation                           | □ Two-stage ■ Immediate placement/Immediate loading                         |
| Methods         |   | □ Dehiscence □ Fenestration □ Horizontal bone loss                      | ■ Intra-socket  □ Bone augmentation  □ Ridge Split            | □ Two-stage ■ Immediate placement/Immediate loading ■ Simultaneous approach |
| Methods         |   | □ Dehiscence □ Fenestration □ Horizontal bone loss □ Vertical bone loss | ■ Intra-socket  □ Bone augmentation  □ Ridge Split  □ Lateral | □ Two-stage ■ Immediate placement/Immediate loading ■ Simultaneous approach |

#### **Before**




The patient visited for discomfort on the maxillary left central and

#### **After**



The periapical X-ray showed a well-integrated implant with a cantilever bridge.

#### Conclusion



and discomfort on the maxillary left lateral incisor.

maxillary left central incisor. Crown fracture of the maxillary left lateral incisor.

3. Pre-operation. CBCT findings.



4. Intra-operation. Extraction and immediate implant placement on the maxillary left central incisor. The implant was positioned on the palatal slope

5. Intra-operation. Extraction of the maxillary left

6. Intra-operation. THE Graft Collagen was used for A healing abutment was connected.



7. Intra-operation. Ridge preservation was done at the 8. Post-operation. Flap approximation with 4-0 nylon. 9. Post-operation, CBCT findings. maxillary left lateral incisor using THE Graft Collagen.







10. Final prosthesis, 3 months and 3 weeks post-operation.

11. Final prosthesis, 3 months and 3 weeks post-operation.

12. Final prosthesis, 3 months and 3 weeks post-operation.

## Gap filling with THE Graft Collagen at the time of immediate implantation

#### Sang-Chul Lim

• Mi-Sa Well Dental Clinic, Korea



#### **Case Summary**

| Nationality     | nd Korean  | Age                        | Late Thirties                 |   |                               |
|-----------------|--|----------------------------|-------------------------------|---|-------------------------------|
| Chief Complaint | Presence of pus in the gu  | ms and a cavity at the too | -                             |   |                               |
| Treatment Plan  | Immediate implant placement after tooth extraction.     Implantation of graft materials into the gap between the placed implant and the extraction socket. |                            |                               |   | 5 4                           |
|                 | Category   | Products                   | Method                        | De  | escription of the method      |
| Materials       | Bone graft   | THE Graft Collagen         | Single used                   | Intra-socket grafting takes place immediately following extraction and implant placement. |                               |
|                 | Membrane   | Not used                   |                               |   |                               |
|                 | Category   | Indication                 | Approach                      |   | Surgical Procedures           |
|                 | ■ Implantology   | ■ Extraction sockets       | ☐ Alveolar ridge preservation | ■ One-stag  | е                             |
|                 | □ Periodontology   | □ Dehiscence               | ■ Intra-socket                | ☐ Two-stage   | е                             |
|                 |  | ☐ Fenestration             | ☐ Bone augmentation           | □ Immediate   | e placement/Immediate loading |
| Methods         |  | ☐ Horizontal bone loss     | □ Ridge Split                 | ■ Simultane   | eous approach                 |
|                 |  | ☐ Vertical bone loss       | □ Lateral                     | □ Staged/D  | elayed approach               |
|                 |  | ☐ Sinus pneumatization     | □ Crestal                     |   |                               |
|                 |  | □ Peri-implantitis         |                               |   |                               |
|                 |  | □ Furcation                |                               |   |                               |



#### **Before**

After

Immediately after extraction.

4 months after intra-socket grafting.

#### Conclusion

In the case of immediate implant placement after extraction, intra-socket grafting could be performed simply by using THE Graft Collagen instead of particle-type bone graft material. Furthermore, by using THE Graft Collagen, particles were fixed at the applied position without scattering, thereby enhancing surgical manipulation and reducing surgery time.

### **Case Summary**



the maxillary first and second premolars.

EXTRACTION SOCKETS | INTRA-SOCKET

- the maxillary first and second premolars.
- and second premolars.



- **4. Intra-operation.** Immediate implant placement at the extraction sockets.
- 5. Intra-operation. After implant placement, healing abutments were connected. And gap filling was performed using THE Graft Collagen.
- 6. Intra-operation. Sutured.



- 7. Post-operation. The day of implant placement.
- 8. Suture removal was performed, 2 weeks after the bone grafting. Good healing was achieved without swelling, pain, or complications.

  9. Post-operation, 4 months. The gingiva around the healing abutments healed well.



- 10. Final prosthesis, 4 months after the bone grafting.

  Clinically observed excellent implant fixation due to

  New bone formation and stable marginal bone levels

  New bone formation and stable marginal bone levels

  Elevated sinus floor and marginal bone levels we could be observed on the panoramic radiograph.
  - Elevated sinus floor and marginal bone levels were



#### HORIZONTAL BONE LOSS | BONE AUGMENTATION

## Contour augmentation at the molar site with early implantation

### Sang-Chul Lim

• Mi-Sa Well Dental Clinic, Korea



#### **Case Summary**

| Nationality     | Korean   | Age                    | Early Sixties                 |   |                               |
|-----------------|--|------------------------|-------------------------------|---|-------------------------------|
| Chief Complaint | Tooth mobility and pain on biting at the lower left area.  |                        |                               |   |                               |
| Treatment Plan  | Extraction of the madibular left first molar due to vertical root fracture.     Re-evaluation after 3-4 months due to periapical pathology.     Implant placement and bone augmentation at the mandibular left first molar region. |                        |                               |   | 6                             |
|                 | Category   | Products               | Method                        | De  | escription of the method      |
| Materials       | Bone graft   | THE Graft Collagen     | Single used                   | In order to achieve buccal contour augmentati<br>THE Graft Collagen was applied on the bucca<br>side. |                               |
|                 | Membrane   | Not used               |                               |   |                               |
|                 | Category   | Indication             | Approach                      |   | Surgical Procedures           |
|                 | ■ Implantology   | ☐ Extraction sockets   | ☐ Alveolar ridge preservation | ■ One-stag  | е                             |
|                 | □ Periodontology   | □ Dehiscence           | □ Intra-socket                | □ Two-stag  | е                             |
|                 |  | □ Fenestration         | ■ Bone augmentation           | □ Immediat  | e placement/Immediate loading |
| Methods         |  | ■ Horizontal bone loss | □ Ridge Split                 | □ Simultane   | eous approach                 |
|                 |  | ☐ Vertical bone loss   | □ Lateral                     | ■ Staged/D  | elayed approach               |
|                 |  | ☐ Sinus pneumatization | □ Crestal                     |   |                               |
|                 |  | ☐ Peri-implantitis     |                               |   |                               |
|                 |  | □ Furcation            |                               |   |                               |



#### Before

Implant placement was performed 3 months after the extraction.



#### **After**

2 months post-GBR. The gingiva around healing abutments healed very well.

#### Conclusion

Buccal contour augmentation can be performed successfully using THE Graft Collagen with stable long-term loutcomes.



1. Post-extraction, healing socket.

Post-extraction, 3 months. Remodeling of the extraction socket

3. Intra-operation. The defect remained partially.



Intra-operation. Implant placement (one-stage) and the healing abutments were connected.

5. Intra-operation. THE Graft Collagen was applied to the buccal side for buccal contour augmentation.

**6. Intra-operation.** THE Graft Collagen was applied and condensed on the buccal side.



7. Intra-operation. Sutured.

8. Post-operation. The implant was ideally positioned. 9. Post-operation, 2 months post-GBR. The gingiva

 Post-operation, 2 months post-GBR. The gingiva around healing abutments healed very well with good volume.



 Final prosthesis, 2 months post-GBR. It was observed that the buccal contour was maintained.

11. Final prosthesis, 2 months post-GBR.
Stable marginal bone.

12. Follow up, 6 months post-GBR. It was observed that the buccal contour was well maintained.

#### HORIZONTAL BONE LOSS | BONE AUGMENTATION

## Contour augmentation for the long-term success of implantation

## Sang-Chul Lim • Mi-Sa Well Dental Clinic, Korea



#### **Case Summary**

| Nationality     | Korean  | Age                                    | Early sixties                 |   |                       |      |
|-----------------|---|--|-------------------------------|---|-----------------------|------|
| Chief Complaint | Mobility in the mandibular anterior teeth.  |  |                               |   |                       |      |
| Treatment Plan  | 1. Extraction of the lower anterior teeth. 2. Implant placement on the mandibular right and left lateral incisors. 3. Layered bone grafting (First layer: autogenous bone, second layer: THE Graft Collagen). 4. Second operation and additional bone graft with THE Graft Collagen. 5. Final prosthesis. |  |                               |   | 2 1                   | 1 2  |
|                 | Category  | Products                               | Method                        | De  | escription of the met | thod |
| Materials       | Bone graft  | Autogenous bone,<br>THE Graft Collagen | Layering                      | The autogenous bone is layered first, followed by THE Graft Collagen. 2 months after the initial bone graft, an additional bone graft was performed to recover from labial bone resorptio due to periodontitis. |                       |      |
|                 | Category  | Indication                             | Approach                      |   | Surgical Procedure    | es   |
|                 | ■ Implantology  | ☐ Extraction sockets                   | ☐ Alveolar ridge preservation | □ One-stage   | е                     |      |
|                 | □ Periodontology  | □ Dehiscence                           | □ Intra-socket                | ■ Two-stage   | е                     |      |
|                 |   | □ Fenestration                         | □ Immediate                   | e placement/Immedia   | ate loading           |      |
| Methods         |   | ■ Horizontal bone loss                 | ☐ Ridge Split                 | ■ Simultane   | eous approach         |      |
|                 |   | □ Vertical bone loss                   | □ Lateral                     | ■ Staged/D  | elayed approach       |      |
|                 |   | ☐ Sinus pneumatization                 | □ Crestal                     |   |                       |      |
|                 |   | □ Peri-implantitis                     |                               |   |                       |      |
|                 |   | □ Furcation                            |                               |   |                       |      |



#### **Before**

Periodontally compromised anterior teeth with gum bleeding in the mandibular anterior area.



#### After

Despite the severe bone loss and subsequent bone defect, the labial contour could be reconstructed by

#### Conclusion

THE Graft Collagen could be successfully used for regeneration of the alveolar ridge. Additionally, the quality of the newly formed bone was sufficient for implant placement.

THE Graft Collagen was easier to manipulate compared to the particulate bone graft material for guided bone regeneration.

#### **Case Summary**



- 1. Pre-extraction. A narrow ridge was observed in the 2. Post-extraction, 1 month. mandibular anterior area.
- 3. Intra-first operation. Insufficient healing of the extraction socket could be seen after flap elevation.



- Intra-first operation. Implant placement on the mandibular right and left lateral incisors.
- **5. Intra-first operation.** Autogenous bone was first layered on the extraction sockets.
- **6. Intra-first operation.** THE Graft Collagen was second-layered on autogenous bone.



7. Intra-first operation.

- 8. Intra-first operation. Flap repositioning and suture. 9. Intra-second operation, 2 months after the bone grafting.
  - The extraction socket cavity was filled with new bone, and the quality of the new bone was very similar to that of the host bone, and there was no visible boundary between them. The depression disappeared but the labial bone was resorbed, and additional bone grafting was performed.



10. Post-first operation, 5 months after the bone grafting. 11. Final prosthesis, 5 months after the bone grafting. 12. Follow up, 7 months after the bone grafting. Volume recovery by creating sufficient bone on the labial side.

The bone level was well maintained

## Space provision in the lower posterior region using THE Graft Collagen

#### Sang-Chul Lim

• Mi-Sa Well Dental Clinic, Korea



#### **Case Summary**

| Nationality     | Korean  | Age  | Mid-Fifties                   |   |
|-----------------|---|--|-------------------------------|---|
| Chief Complaint | Difficulty in chewing becar                         | use of the mobility of low   |                               |   |
| Treatment Plan  | due to failed prostheses 2. Implantation of the man | Extraction of the mandibular left first and second molars due to failed prostheses.     Implantation of the mandibular left first and second molars with GBR on the buccal site for horizontal augmentation. | Grafting Area 6 7             |   |
|                 | Category  | Products   | Method                        | Description of the method   |
| Materials       | Bone graft  | THE Graft Collagen   | Single used                   | GBR was performed to increase the narrowed ridge on the buccal-lingual side. THE Graft Collagen was used to make space. |
|                 | Membrane  | Collagen membrane  | Primary closing               | Buccal bone defect covered.   |
|                 | Category  | Indication   | Approach                      | Surgical Procedures   |
|                 | ■ Implantology                                      | ☐ Extraction sockets   | ☐ Alveolar ridge preservation | ■ One-stage   |
|                 | □ Periodontology                                    | □ Dehiscence   | □ Intra-socket                | ☐ Two-stage   |
|                 |   | □ Fenestration   | ■ Bone augmentation           | ☐ Immediate placement/Immediate loading   |
| Methods         |   | ■ Horizontal bone loss   | ☐ Ridge Split                 | ■ Simultaneous approach   |
|                 |   | ■ Vertical bone loss   | □ Lateral                     | ■ Staged/Delayed approach   |
|                 |   | ☐ Sinus pneumatization   | □ Crestal                     |   |
|                 |   | □ Peri-implantitis   |                               |   |
|                 |   | ☐ Furcation  |                               |   |



#### **Before**

4 months after extraction, implant

After

2 months post-GBR, final prosthesis.

#### Conclusion

When GBR is performed at the same time as implant placement, it is thought that THE Graft Collagen is more advantageous in making the space than the particulate bone. Furthermore, THE Graft-THE Graft Collagen combination technique was thought to be a simple space-creating device.



1. Pre-extraction.

2. Post-extraction, 5 weeks.

3. Post-extraction, 4 months.



Intra-operation. Implant placement (one-stage) and healing abutments were connected.

Intra-operation. Simultaneously with implant placement, GBR was performed to increase the narrowed ridge on the buccal-lingual side. THE Graft Collagen was applied to maintain space.

**6. Intra-operation.** Covered with a resorbable collagen membrane.



7. Intra-operation. Sutured.

8. Post-operation. Well-positioned implants were

9. Final prosthesis, 2 months post-GBR.



10. Final prosthesis, 2 months post-GBR.

11. Final prosthesis, 2 months post-GBR.

Follow up, 5 months post-GBR.
 Marginal bone and crestal bones are remodeled and stable. The buccal ridge remained stable, so the soft tissues also remained stable.

## Minimalistic bone augmentation procedure using a combination of THE Graft and THE Graft Collagen


#### Sang-Chul Lim

• Mi-Sa Well Dental Clinic, Korea



#### **Case Summary**

| Nationality     | Korean   | Age                              | Early forties                 |  |                               |
|-----------------|--|----------------------------------|-------------------------------|--|-------------------------------|
| Chief Complaint | Treatment for the front teeth was required.  |                                  |                               |  |                               |
| Treatment Plan  | 1. Removal and evaluation of the cantilever bridge on the maxillary left and right central incisors. 2. Implantation and bone grating were done on the maxillary left lateral incisor pontic area simultaneously. 3. Immediate implantation with gap filling was performed on the maxillary left and right central incisors. 4. Final restoration is supported by the maxillary right central incisor, left central and lateral incisors implants. |                                  |                               | Grafting<br>Area   | 2                             |
|                 | Category   | Products                         | Method                        | De   | escription of the method      |
| Materials       | Bone graft   | THE Graft,<br>THE Graft Collagen | Layering                      | To restore the narrow ridge, THE Graft and THE Graft Collagen were applied using the layering technique. |                               |
|                 | Category   | Indication                       | Approach                      |  | Surgical Procedures           |
|                 | ■ Implantology   | ☐ Extraction sockets             | ☐ Alveolar ridge preservation | □ One-stag   | е                             |
|                 | □ Periodontology   | □ Dehiscence                     | □ Intra-socket                | ■ Two-stag   | е                             |
|                 |  | □ Fenestration                   | ■ Bone augmentation           | □ Immediat   | e placement/Immediate loading |
| Methods         |  | ■ Horizontal bone loss           | □ Ridge Split                 | ■ Simultane  | eous approach                 |
|                 |  | ■ Vertical bone loss             | □ Lateral                     | ■ Staged/D   | elayed approach               |
|                 |  | ☐ Sinus pneumatization           | □ Crestal                     |  |                               |
|                 |  | □ Peri-implantitis               |                               |  |                               |
|                 |  | □ Furcation                      |                               |  |                               |



#### **Before**

Observation of the narrow ridge of the maxillary left lateral incisor.



#### After

A new bone was formed very well up to the top of the cover screw, and the labial volume was recovered. It was the new bone and the host bone.

#### Conclusion


When THE Graft Collagen was applied, the sufficient new bone formation could be confirmed even in narrow ridges. Bone regeneration went well up to the top of the cover screw. The quality of the bone was good enough and did not loosen during drilling. The healing was very good, and after the new bone was formed, the particles of the bone graft material did not move in the soft



1. Pre-operation. The old bridge was moving.

of the maxillary left lateral incisor site (pontic).

2. Intra-first operation. Observation of the narrow ridge 3. Intra-first operation. Observation of the narrow ridge of the maxillary left lateral incisor pontic area.



4. Intra-first operation. Implant placement.

5. Intra-first operation. THE Graft was first layered on the labial area.

**6. Intra-first operation.** THE Graft Collagen was layered on top of THE Graft.



7. Intra-first operation. Sutured.

8. Post-first operation. A temporary bridge was built. 9. Intra-second operation, 5 months after the bone grafting.

A new bone was formed very well up to the top of the cover screw, and the labial volume was recovered. It was difficult to check the boundary between the new bone and the host bone.



10. Final prosthesis, 9 months after the bone grafting. 11. Final prosthesis, 9 months after the bone grafting. 12. Follow up, 1 year and 8 months after the bone grafting.

## Overcoming sinus membrane perforation with THE Graft Collagen

#### Sang-Chul Lim

• Mi-Sa Well Dental Clinic, Korea



#### **Case Summary**

| Nationality     | Korean  | Age                    | Late fifties                  |  |  |  |
|-----------------|---|------------------------|-------------------------------|--|--|--|
| Chief Complaint | Difficulty in chewing because of the absence of the maxillary right first premolar and first molar. |                        |                               |  |  |  |
| Treatment Plan  | Sinus elevation with imp     GBR on narrow ridge si     Zirconia bridge.                            |                        | Grafting<br>Area              | 6 4  |  |  |
| Materials       | Category  | Products               | Method                        | Description of the method  |  |  |
|                 | Bone graft  | THE Graft Collagen     | Single used                   | THE Graft Collagen was applied without a membrane covering the maxillary sinus perforation. Following that, GBR was performed with THE Graft Collagen. |  |  |
|                 | Membrane  | Not used               |                               |  |  |  |
| Methods         | Category  | Indication             | Approach                      | Surgical Procedures  |  |  |
|                 | ■ Implantology  | ☐ Extraction sockets   | ☐ Alveolar ridge preservation | □ One-stage  |  |  |
|                 | □ Periodontology  | □ Dehiscence           | □ Intra-socket                | ■ Two-stage  |  |  |
|                 |   | ☐ Fenestration         | ■ Bone augmentation           | ☐ Immediate placement/Immediate loading  |  |  |
|                 |   | ☐ Horizontal bone loss | □ Ridge Split                 | ■ Simultaneous approach  |  |  |
|                 |   | ☐ Vertical bone loss   | ■ Lateral                     | ☐ Staged/Delayed approach  |  |  |
|                 |   | ■ Sinus pneumatization | □ Crestal                     |  |  |  |
|                 |   | □ Peri-implantitis     |                               |  |  |  |
|                 |   | □ Furcation            |                               |  |  |  |



#### **Before**




#### After

5 months after sinus augmentation, final

#### Conclusion

During maxillary sinus lifting, the maxillary sinus membrane was often torn. At this time, it was better to use THE Graft Collagen than particulate bone graft material. THE Graft Collagen did not flow through the perforated sinus membrane, which was the cause of the

It may be suggested that without covering the perforated sinus membrane with a membrane, it was possible to perform a sufficiently successful operation with THE Graft Collagen alone.



1. Pre-operation

2. Pre-operation.


3. Intra-first operation. Maxillary sinus membrane perforation during maxillary sinus detaching.



4. Intra-first operation. THE Graft Collagen, which was less likely to scatter than particulate bone material

5. Intra-first operation. THE Graft Collagen was used for simultaneous GBR and implant placement.

6. Intra-first operation. Covered with a resorbable simultaneous GBR and implant placement. was applied for sinus augmentation.



7. Intra-first operation. Sutured.

8. Intra-first operation. Well-placed implants were observed, and THE Graft Collagen applied in the maxillary sinus was observed.

9. Intra-second operation, 5 months post-GBR.



10. Intra-second operation, 5 months post-GBR. The healing abutments were connected and sutured.

Final prosthesis, 6 months post-GBR.
 Delivered zirconia bridge.

12. Final prosthesis, 6 months post-GBR. New bone formation was observed in the maxillary sinus.

#### PERI-IMPLANTITIS | BONE AUGMENTATION

## Regenerative therapy of peri-implantitis in the mandibular first molar area

#### Yong-Seok Cho

• The 22<sup>nd</sup> Century Seoul Dental Hospital



#### **Case Summary**

| Nationality     | Korean                    | Age   | Early Sixties                 |   |  |  |
|-----------------|---------------------------|---|-------------------------------|---|--|--|
| Chief Complaint | Gingival swelling and ble | eding around the mandibul                                 | Grafting                      |   |  |  |
| Treatment Plan  |                           | e mandibular left first mola<br>d bone and non-resorbable | Area                          | 6   |  |  |
| Materials       | Category                  | Products  | Method                        | Description of the method   |  |  |
|                 | Bone graft                | THE Graft Collagen  | Single used                   | GBR with collagenated bone and non-resorbable membrane.   |  |  |
|                 | Membrane                  | OpenTex-TR  | Primary closure               | Covered with a non-resorbable titanium-reinforced polytetra-fluorethylene (PTFE) membrane over the graft. |  |  |
|                 | Others                    | HybenX  |                               | Chemical cleaning of the mandibular left first molar implant surface using HybenX.                        |  |  |
| Methods         | Category                  | Indication  | Approach                      | Surgical Procedures   |  |  |
|                 | ■ Implantology            | □ Extraction sockets                                      | ☐ Alveolar ridge preservation | □ One-stage   |  |  |
|                 | □ Periodontology          | □ Dehiscence  | □ Intra-socket                | ☐ Two-stage   |  |  |
|                 |                           | □ Fenestration  | ■ Bone augmentation           | ☐ Immediate placement/Immediate loading   |  |  |
|                 |                           | ☐ Horizontal bone loss                                    | □ Ridge Split                 | ☐ Simultaneous approach   |  |  |
|                 |                           | □ Vertical bone loss                                      | □ Lateral                     | ☐ Staged/Delayed approach   |  |  |
|                 |                           | ☐ Sinus pneumatization                                    | □ Crestal                     | ☐ Staged/Delayed approac  |  |  |
|                 |                           | ■ Peri-implantitis  |                               | Peri-implantitis treatment by mechanical an chemical cleaning.  |  |  |
|                 |                           | □ Furcation   |                               |   |  |  |



#### **Before**

The patient visited for discomfort and bone loss around the mandibular left first molar implant.



#### After

A periapical X-ray showed well-regenerated bone around the mandibular left first molar implant.

#### Conclusion

GBR using THE Graft Collagen and a non-resorbable membrane resulted in successful bone regeneration in the implant with peri-implantitis.







around the implant of the mandibular left first molar, which indicated peri-implantitis.




3. Intra-first operation. Mechanical cleaning of the mandibular left first molar implant surface using the scaler. Chemical cleaning of the mandibular left first molar implant surface using HybenX.





4. Intra-first operation. THE Graft Collagen was grafted around the mandibular left first molar implant.

5. Intra-first operation. Over the graft, a non-resorbable 6. Intra-first operation. A flap approximation was done. titanium-reinforced PTFE membrane was covered with OpenTex-TR.



7. Post-first operation.



8. Post-first operation, 1 month post-GBR.



9. Post-second operation, 2 months post-GBR.



Post-operation, 3 months and 2 weeks post-GBR. Gingival attachment failure around implants. So, a free gingival graft was done.



11. Post-operation, 1 year and 9 months post-GBR. A healthy gingival condition was observed.



**12. Post-operation, 1 year and 9 months b. post-GBR.** A CT view showed bone regeneration around the mandibular left first molar implant.

## Author

## Prof. Jae-Hong Lee

Wonkwang University

DDS, MSD, Ph.D, Yonsei University College of Dentistry,
Periodontist, Diplomate Korean Board of Periodontology
Clinical and Research Fellow, Yonsei University Dental Hospital
Clinical and Research Fellow, NHIS Ilasn Hospital
Associate Professor, Department of Periodontology, Daejeon Dental Hospital,
Wonkwang University College of Dentistry



## Prof. Jung-Seok Lee

Yonsei University

DDS, MSD, Ph.D, Yonsei University College of Dentistry
Periodontist, Diplomate Korean Board of Periodontology
Chair, Department of Periodontology, Yonsei University College of Dentistry
Director, Department of Periodontology, Yonsei University Dental Hospital
Visiting Scholar, Department of Oral Biology, Medical University of Vienna, Austria



## Prof. Dong-Woon Lee

Wonkwang University

DDS, MSD, Ph.D, Wonkwang University,
Periodontist, Diplomate Korean Board of Periodontology
Associate Professor, Department of Periodontology, Wonkwang University
College of Dentistry, Iksan, Korea
Residency, Department of Periodontology, Veterans Health Service Medical

Center, Seoul, Korea
Chair and Faculty, Department of Periodontology, Veterans Health Service

Medical Center, Seoul, Korea

Visiting Scholar, Department of Advanced periodontics, University of Southern California, Los Angeles, USA

Clinical Associate Professor, Department of Periodontology, Dental Hospital, Wonkwang University College of Dentistry, Iksan, Korea



## Dr. Yong-Seok Cho

22nd Century Seoul Dental Hospital

DDS, Ph.D, Seoul National University,
Diplomate, Korean Board of Oral & Maxillofacial Surgery and Korean Board of
Maxillofacial Plastic and Reconstructive Surgery
Chairman, Department of Dentistry, Chungbuk National University Hospital
Adjunct professor, Seoul National University School of Dentistry
Adjunct professor, Asan Medical Centre College
Adjunct professor, Medicine and Ulsan University
Director, AIC (Apsun Dental Implant Research and Education Centre)



## Dr. Sang-Chul Lim

Director, The 22nd Century Seoul Dental Hospital

WELL Dental Clinic at Misha

DMD, Ph.D, Kyung Hee University,
Periodontist, Diplomate Korean Board of Periodontology
Visiting Professor, University of Toronto, Faculty of Dentistry, Toronto, Canada
Visiting Scientist, University of British Columbia, Faculty of Dentistry, Vancouver,
Canada

Adjunct Professor, Kyung Hee University Dental Hospital Director, WELL Dental Clinic



Editor-in-Chief Prof. Yeek Herr

Editorial Boar

Dr. Kyoung-Man Min, Dr. Chang-Kyun Lee, Prof. Jin-Young Park